




CHARACTERIZING SELF ASSEMBLED MONOLAYERS USING SCANNING TUNNELING MICROSCOPY

By Robert Conwell Advisor: Dr. Lloyd Bumm 30 July 2019

SCANNING TUNNELING MICROSCOPE (STM) FUNDAMENTALS


- Electrons tunnel from atomically sharp tip to flat conductive sample $I \propto e^{-kd}$
- Set point current and negative feedback for constant height
 - Low I: Tip too far, move forward
 - High *I*:Tip too close, move back
- Images surface topography for homogeneous samples


Piezo Actuators





STM IMAGING EFFECTS

- Various imaging effects which distort the image
 - Drift
 - Hysteresis
 - Creep
- Either image is distorted or SAMs are not a lattice
 - Latter is highly unlikely
- Apply image correction in MATLAB











M.P.Yothers, A. E. Browder, and L.A. Bumm, Rev. Sci. Instrum. 88 (1), 013708 (2017).




M.P.Yothers, A. E. Browder, and L.A. Bumm, Rev. Sci. Instrum. 88 (1), 013708 (2017).



M.P.Yothers, A. E. Browder, and L.A. Bumm, Rev. Sci. Instrum. 88 (1), 013708 (2017).

ALKANETHIOL SELF ASSEMBLED MONOLAYERS (SAMs)

- Alkanethiol hydrocarbon chain with sulfur at the end
- Self assemble spontaneously grow 2D crystal from alkanethiol vapor
- Monolayer Only a molecule thick




Decanethiol – CI0

ALKANETHIOL SELF ASSEMBLED MONOLAYERS (SAMs)

- Alkanethiol SAMs on gold substrate
- Abbreviate name by length of hydrocarbon chain

Name	Abbreviation
Octanethiol	C8
Nonanethiol	С9
Decanethiol	C10
Undecanethiol	CII
Dodecanethiol	C12



C10 Sample 2019 #3 2019_06_17_0033


CRYSTAL LATTICE

- Crystal repeating arrangement of atoms
- Two components
 - Lattice \rightarrow How to repeat
 - Basis \rightarrow What is being repeated




 Unit cell → Unit cell is the smallest group of atoms which can be repeated to reconstruct the entire crystal, i.e. contains all information of the crystal

UNIT CELL & LOCATION CONFIDENCE



- Averaged unit cell image
 - Grayscale
- Ellipses location confidence
 - Colored by basis site
 - Center Mean location
 - Inner ellipse $I\sigma$, outer 2σ


UNIT CELL & LOCATION CONFIDENCE




- Averaged unit cell image
 - Grayscale
- Ellipses location confidence
 - Colored by basis site
 - Center Mean location
 - Inner ellipse $I\sigma$, outer 2σ

ALKANETHIOL CRYSTAL


- All 4 molecule colors are alkanethiol
- Structure is very similar to a one-molecule basis
- How do we compare these structures?




- Molecule ends are close to one-molecule basis lattice sites
- Fingerprint plot shows distance from these sites
- Easier to see differences between structures




- Molecule ends are close to one-molecule basis lattice sites
- Fingerprint plot shows distance from these sites
- Easier to see differences between structures





- Molecule ends are close to one-molecule basis lattice sites
- Fingerprint plot shows distance from these sites
- Easier to see differences between structures





- Molecule ends are close to one-molecule basis lattice sites
- Fingerprint plot shows distance from these sites
- Easier to see differences between structures





MOTIVATION

Does selecting molecules on/near a boundary skew results?

- Close to a boundary, the crystal may be strained
- Analysis on boundary vs Analysis 5-10 molecules away from boundary


- Are chains of odd and even length structurally different?
 - Can we tell odd from even by looking at them?


MOTIVATION

Does selecting molecules on/near a boundary skew results?


- Close to a boundary, the crystal may be strained
- Analysis on boundary vs Analysis 5-10 molecules away from boundary

Are chains of odd and even length structurally different?
Can we tell odd from even by looking at them?







Normal Analysis




On-boundary Analysis





Offset plots





2019_07_17_0069 C8/C9 (2019 #6)






2019_07_17_0069 C8/C9 (2019 #6)

0.04 0.03 0.02 0.01 Distance (nm) 0 -0.01 -0.02 0 -0.03 -0.04 -0.05 -0.04 -0.03 -0.02 -0.01 0 0.01 0.02 0.03 0.04 Distance (nm)

Normal Analysis

On-boundary Analysis



Offset plots

BOUNDARY ANALYSIS

• Confidence ellipses are larger

Increases the uncertainty

• Exclude molecules near the edge


MOTIVATION


• Does selecting molecules on/near a boundary skew results?


- Close to a boundary, the crystal may be strained
- Analysis on boundary vs Analysis 5-10 molecules away from boundary

• Are chains of odd and even length structurally different?


• Can we tell odd from even by looking at them?










STRUCTURE COMPARISON – CI0









Distance (nm)

CONCLUSIONS AND FUTURE DIRECTION

- Made Alkanethiol SAM samples and imaged with STM
- Measured effects of proximity when analyzing images with boundaries and found noticeable effect
- Compared fingerprints of odd and even chain length SAMs and did not find a difference

ACKNOWLEDGMENTS

- Dr. Lloyd Bumm, Mitch Yothers, and Bumm group
- Dr. Abbott, Dr. Strauss, Amber Roepe, organizers
- National Science Foundation
- Nvidia